



## Minimum System Requirements \*

### LABTrack Personal and Personal Lite Editions

Windows XP Professional, Vista Professional, Business or Ultimate (Home versions of Windows are not supported)  
Pentium 4, 1 GHZ CPU or better  
2 GB RAM, 20 GB free disk space  
.NET Framework 3.5  
IIS 5.0 or newer  
Internet Explorer 6, 7 or 8, Firefox 2 or 3 or Safari 3.2.  
Either Oracle XE or SQL Express (provided with LABTrack)

### LABTrack Pro Team Edition

Windows 2000 Server OS or newer  
Pentium 4, 2.0 GHZ CPU or better  
2 GB RAM or more, 20 GB free disk space  
Either SQL Server 2005 or newer or Oracle 10g  
IIS 5.0 or newer  
.NET Framework 3.5  
Internet Explorer 6, 7 or 8, Firefox 2 or 3 or Safari 3.2.  
Microsoft Reporting Services or Crystal Reports XI Web Server  
Version required for custom reports.  
Application and Database Server can be hosted on the same server.  
Additional disk space as required by database storage.

### LABTrack Pro Enterprise Edition

**Application Server** - Windows 2000 Server OS or newer  
Pentium 4, 2.0 GHZ CPU or Better, 10 GB free disk space  
4 GB RAM or more  
IIS 5.0 or newer  
.NET Framework 3.5  
Microsoft Reporting Services or Crystal Reports XI Web Server  
Version required for custom reports.  
Internet Explorer 6 or 7, Firefox 2 or 3 or Safari 3.2.

**Database Server** – Windows 2000 Server OS or newer  
and either SQL Server 2005 or Oracle 10g

- or -

Linux, AIX or Solaris OS and Oracle 10g  
Additional disk space as required by database storage.

### Team & Enterprise Client PC Requirements

Windows XP, Vista or Mac X OS or newer  
1 GB RAM, 50 MB free disk space  
Internet Explorer 6, 7 or 8, Firefox 2 or 3 or Safari 3.2.

\* Requirements are subject to change without notice.

**Electronic Lab Notebook (ELN)** -  
Designed to replace paper notebooks for Research, Analytical and Production laboratories. LABTrack meets all of the legal, regulatory and typical corporate requirements for maintaining laboratory notebooks.

The **LABTrack Pro Team Edition** is designed and priced attractively for smaller groups numbering from 5 to 200 users.

The **LABTrack Pro Enterprise Edition** is designed for larger scale implementations numbering from 50 to 10,000 users. The Enterprise Version is developed to support typical WAN bandwidth and global requirements.

The **LABTrack Personal Edition** installs on a Windows PC allowing users to edit notebooks locally. When combined with a Team or Enterprise version, users may check-in and check-out notebooks from the master database.

---

## LABTrack, LLC

22365 El Toro Rd #185, Lake Forest, CA 92630  
Tele: (US) +1 (949) 528-3840, [www.labtrack.com](http://www.labtrack.com)

