

Introducing LABTrack Personal/Q1 the Electronic Lab Notebook Appliance

Introductory Price - Only \$ 2,500

LABTrack Personal/Q1 can bring your laboratory notebook into the 21st century. Incorporating all the features of the award winning Team edition of the LABTrack Software – LABTrack Personal/Q1 ELN provides all the flexibility of paper and the portability of an ultra mobile PC with the searching and file capture capability of this amazing electronic lab notebook application. Up to 5 users can share this Personal ELN with unique user accounts.

As a stand alone appliance, LABTrack Personal/Q1 gives scientists and engineers the ability to replace their paper notebooks with a powerful electronic equivalent. Like the Team and Enterprise editions, LABTrack Personal/Q1 provides an easy to learn word processor coupled with an advanced Relational Database. Up to 2 GB of notebook data can be stored in the included database.

When combined with either the Team or Enterprise editions of LABTrack, Personal/Q1 users can check-out and check-in notebooks from a master LABTrack server database.

The Personal/Q1's stylus and built-in software keyboard makes data entry very easy. Add a monitor, a USB keyboard and a mouse and your Q1 becomes a fully functional Windows desktop PC.

It has never been easier to capture and store all of your scientific data. LABTrack Personal/Q1 is where *Science is stored in the 21st century™*.



Product Features and Technical Details

- **Pre-Installed LABTrack Personal ELN Software**
- **Pre-Installed MSDE Database** – 2 GB limit.
- **Samsung Q1 Ultra Mobile PC**
- **Full-powered PC** with multimedia functionality
- **Weighs just 1.7 pounds**
- **7-inch LCD touchscreen** (800 x 480 pixels) with landscape and portrait orientation
- **40 GB hard disk drive; 512 MB of RAM** (1 GB max)

- **Integrated WLAN 802.11 b/g and Bluetooth 2.0** to keep you connected anytime, anywhere.

- Powered by **Windows XP Tablet and 1.0 GHz Intel Pentium M ULV processor**

Q1 Processor, Memory, and Motherboard

- **Hardware Platform:** PC
- **Processor:** 900 MHz
- **System Bus Speed:** 400
- **RAM:** 512 MB
- **RAM Type:** DDR2 SDRAM

Dimensions

- **Size (LWH):** 11.3 inches, 7.8 inches, 4.13 inches

Power

- **Rated Charge (normal use):** 3 hours

Warranty and Support

- **Includes 90 days telephone and email support.**

Ordering Information

Fax purchase orders to 208-498-3840 attn: P/Q1 Introductory Offer

Order On-Line with Major Credit Cards by browsing to www.Labtrack.com

Discounts Available on Multiple Unit Orders. Please call 949-528-3840 for pricing.



LABTrack, LLC | 22365 El Toro Rd #185 | Lake Forest, CA 92630 | Tele: 949-528-3840
Fax: 208-498-7627

Copyright © 2007, LABTrack LLC. All Rights Reserved. LABTrack and LABTrack Personal are trademarks of LABTrack, LLC.