



LABTrack Tech Note 08-011 Personal Edition Admin Menu Fix

Applicable only to LABTrack Personal and Personal Lite Edition rev. 3.0.7 / SQL Express

Because of an incorrect setting in the shipped database, commonly used Administration Wizard and Administration Option menu items may not be visible and activated. In particular the functions for adding new notebooks, users, groups and updating the same are affected.

The following steps correct this configuration issue:

1. Start the Microsoft SQL Server Management Studio Express
2. Log in with the SA account. The default password is 'living' without the quotes.
3. Expand the databases and select the LabtrackPersonal database.
4. Click on the **New Query** toolbar button
5. Select the LabtrackPersonal database from the list to the immediate left of the '! Execute' toolbar button.
6. Copy the following command and paste it into the blank Query in the right hand panel:

```
Update Dbsetup SET LocalInstallation =  
'DA02A00B01577EAE19E7A55126AC9CE0F5741EDB', NotebookCheckoutAvail  
able = 'N';
```

7. Click on the **! Execute** button. The message '1 row updated' should appear if the update was successful.
8. Close SQL Server Management Studio Express.
9. Login to LABTrack Personal / Personal Lite and confirm that the Admin Wizard and Admin Option menus are fully populated with the missing menu items.